


# ナフタレンd8標準物質

GC-MS 用 内部標準物質


SPEX 社は ISO17025 試験所認定を取得し、  
ISO Guide34 に準拠した認証標準物質を販売しております

SPEX CertiPrep.®

Components	CAS#	Concentration
Naphthalene-d8	1146-65-2	1000 µg/mL



**SPEXcertifcate**®  
Certificate of Reference Material



**Catalog Number:** S-2660      **Lot No.:** E1140818013

**Description:** Naphthalene-d8      **Ship Date:** 9-14-2015

**Matrix:** Methanol (Purge & Trap Grade)      **Expiration Date:** 9-13-2018

This SPEXOrganics® Certified Reference Material, CRM, is intended primarily for use as a calibration standard or quality control standard for organic chromatography instrumentation such as GC, GC-MS, LC, and LC-MS. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

Compound	CAS #	Labeled	Purity	Certified†	Uncertainty
Naphthalene-d8	1146-65-2	1000 µg/mL	99.3%	1000 µg/mL	± 35 µg/mL

**Final Solution Verification:**  
Final solution integrity verified by Gas Chromatography/Mass Spectrometry. The mass spectrum of each compound was confirmed against the NIST mass spectral database.

† Certified concentration based on gravimetric weights and corrected for the purity of the compound(s) used to prepare the standard. Analytical balance calibration is verified daily with C1 weight set #23-190006 which is registered with Atlantic Scale, and traceable to NIST and NJ Division of Weights and Measures.

This CRM is guaranteed stable and accurate to within the uncertainty listed for the certified value. This includes uncertainty components due to preparation, homogeneity, short term and long term stability. During the stated period of validity, the purchaser will be notified if this product is recalled due to any significant changes in the stability of the solution. For further information, contact the Sales Support Department at crmsales@spexcorp.com.

**Date of Certification:** 9-14-2015      **Certifying Officer:** *Julian Burton*

Page 1 of 2



## ◆ 製品詳細 ◆

製品コード : S-2660  
濃度 : 1000µg/ml  
容量 : 1ml  
Matrix : メタノール

製品コード : S-2660-ACN  
濃度 : 1000µg/ml  
容量 : 1ml  
Matrix : アセトニトリル

西進商事株式会社

神戸本社      TEL. 078-303-3810(代)  
東京支店      TEL. 03-3459-7491(代)  
名古屋営業所      TEL. 052-586-4741(代)  
北海道営業所      TEL. 011-221-2171(代)